



PROCLEAR TORIC XR: NOW WITH EVEN MORE PARAMETER OPTIONS.

PROCLEAR TORIC: OFFER MORE COMFORT TO EVEN MORE OF YOUR ASTIGMATIC PATIENTS.

ProcLEAR Toric has a proven track record of providing outstanding health and comfort. With two base curves, and now a further expanded XR range, even more of your astigmatic patients can enjoy true all-day comfort.






WITH MORE THAN 45,000 AVAILABLE PARAMETERS, MORE OF YOUR ASTIGMATIC PATIENTS CAN ENJOY PROCLEAR'S ALL-DAY COMFORT.

ALL-DAY COMFORT.

- Proclear's patented PC Technology™ creates a unique lens material that contains phosphorylcholine (PC)
- PC molecules attract and surround themselves with water, keeping the lenses moist, healthy, and comfortable throughout the day
- Features the only lens material with an FDA-cleared indication relating to the improvement of comfort for those who experience dryness or mild discomfort during lens wear¹
- Soft, low modulus offers better comfort and fit



EXCELLENT HEALTH.

- PC Hydrogel™ lenses significantly exceed the minimum level of oxygen for safe daily wear²
- The only monthly toric lens that uses PC Technology to create a unique biocompatible lens material and a deposit-resistant surface that remains clean and healthy

LOWER MODULUS.

The lens material of PC Hydrogel lenses has a significantly lower modulus (less stiff) than silicone hydrogel lenses. And a lower modulus can reduce the risk of adverse effects on the cornea such as SEALs and conjunctival flaps.

A WIDER RANGE OF PARAMETERS.

More than 45,000 available parameters enable you to fit more astigmatic patients more precisely.

	MATERIAL	WEARING SCHEDULE	WATER CONTENT	BASE CURVE (mm)	DIAMETER (mm)	SPHERE POWER	CYLINDER POWER	MFG. PROCESS	AXIS
PROCLEAR TORIC	omafilcon A	Daily Wear	62%	8.4 8.8	14.4	+6.00D to -8.00D (0.50D steps above -6.50D)	-0.75, -1.25, -1.75, -2.25	Cast Molded	10° to 180° (10° steps)
PROCLEAR TORIC XR*	omafilcon A	Daily Wear	59%	8.4 8.8	14.4	+10.00D to -10.00D (0.50D steps above +/- 6.50D)	-2.75 to -5.75 (0.50 steps)	FIPS™	5° to 180° (5° steps)
						+6.25D to +10.00D (0.50D steps after +6.50D) -8.50D to -10.00D (0.50D steps)	-0.75, -1.25, -1.75, -2.25		

*Proclear Toric XR is a made-to-order lens and will ship within 10 days from the time you order.

FOR MORE INFORMATION ABOUT PROCLEAR TORIC, CONTACT YOUR COOPERVISION AREA SALES MANAGER OR COOPERVISION AUTHORIZED DISTRIBUTOR.

coopervision.com

800.341.2020



SEE BEYOND THE ORDINARY®

[†]Some restrictions may apply. See coopervision.com for details.
¹ FDA-cleared labeling indication: "May provide improved comfort for contact lens wearers who experience mild discomfort or symptoms relating to dryness during lens wear." ² Holden BA, Mertz GW. Critical oxygen levels to avoid corneal edema for daily and extended wear contact lenses. Invest Ophthalmol. Vis Sci. 1984;25:1161-7; Harvitt DM, Bonnano JA. Re-evaluation of the oxygen diffusion model for predicting minimum contact lens Dk/t values needed to avoid cornea anoxia. Optom Vis Sci. 1999;76:712-9.
 ©2009 CooperVision, Inc. MF 1/09 COO-PRC-5194