



PRODUCT REFERENCE GUIDE

Look to the people, products, and services of CooperVision.

CONTENTS

Spheres 1 – 7

Torics 8 – 10













Multifocals 11

Customer Service: 1-800-341-2020
www.CooperVision.com







PC
TECHNOLOGYVISIBILITY
TINTLOW
MODULUSNATURALLY
WETTABLEHIGH OXYGEN
TRANSMISSIBILITYASPHERIC
OPTICSOPTIMIZED
BALLAST DESIGNUV
BLOCKING

Sphere Lenses





| BRAND | MATERIAL/ H ₂ O CONTENT | BASE CURVE (mm) | DIA (mm) | SPHERE POWER | FDA Group | WEARING SCHEDULE | Dk/t (at -3.00D) | REVENUE CARTON SIZE | TRIALS LABELED AS/ PACK SIZE | FEATURES |
|--|---------------------------------------|-----------------------|-------------|--|--------------|---------------------|------------------------|---------------------------|---|---|
| 1 Day Disposable Lenses | | | | | | | | | | |
| ClearSight™ 1-Day  | ocufilcon B/ 52% | 8.7 | 14.2 | +6.00D to -10.00D (0.50 steps after -6.00D and +5.00D) | 4 | Daily Disposable | 24 | 90-day and 30-day | 52% 1-Day Trial lenses (ocufilcon B)/ 30 |  |
| Proclear® 1-Day  | omafilcon A/ 60% | 8.7 | 14.2 | +8.00D to -12.00D (0.50 steps after -6.00D and +5.00D) | 2 | Daily Disposable | 36.6 | 90-day and 30-day | 60% 1-Day Trial lenses (omafilcon A)/ 30 |  |
| 2- Week Replacement Lenses | | | | | | | | | | |
| Avaira®  | enfilcon A/ 46% | 8.4 | 14.2 | +0.25 to +8.00D (0.50 steps after +6.00D) | SiHy (1) | Daily | 125 | 6-pack blisters | enfilcon Dk100 Aspheric/ 3 |  |
| | | 8.5 | 14.2 | -0.25D to -12.00D (0.50 steps after -6.00D) | | | | | | |
| Biomedics® 55 UV Premier  | ocufilcon D/ 55% | 8.6 | 14.2 | -0.25D to -10.00D (0.50 steps after -6.00D) | 4 | Extended | 28.1 | 6-pack blisters | 55 UV aspheric contact lens/ 6 |  |
| | | 8.9 | | +0.25D to +6.00D (0.50 steps after +5.00D) | | | | | | |
| Biomedics® 55 UV  | ocufilcon D/ 55% | 8.6 | 14.2 | -0.25D to -10.00D (0.50 steps after -6.00D) | 4 | Extended | 28.1 | 6-pack blisters | 55 UV soft contact lens/ 6 |  |
| | | 8.9 | | +0.25D to +6.00D (0.50 steps after +5.00D) | | | | | | |
| Biomedics® 38  | polymacon/ 38% | 8.6 | 14.0 | -0.25D to -10.00D (0.50 steps after -6.00D) | 1 | Extended | 21 | 6-pack blisters | 38 (polymacon) /6 |  |

PC
TECHNOLOGYVISIBILITY
TINTLOW
MODULUSNATURALLY
WETTABLEHIGH OXYGEN
TRANSMISSIBILITYASPHERIC
OPTICSOPTIMIZED
BALLAST DESIGNUV
BLOCKING

Sphere Lenses

| BRAND | MATERIAL/ H ₂ O CONTENT | BASE CURVE (mm) | DIA (mm) | SPHERE POWER | FDA Group | WEARING SCHEDULE | Dk/t (at -3.00D) | REVENUE CARTON SIZE | TRIALS LABELED AS/ PACK SIZE | FEATURES |
|--|---------------------------------------|-----------------------|-------------|---|--------------|--------------------------------------|------------------------|---------------------------|------------------------------------|---|
| Biomedics® XC  | omafilcon A/ 60% | 8.5 | 14.2 | Plano to -10.00D +0.50D to +6.00D (0.50 steps after -6.00D) | 2 | Daily | 44 | 6-pack blisters | 60%XC (omafilcon A)/ 3 |  |
| Vertex Sphere®  | methafilcon A/ 55% | 8.6 | 14.2 | -0.25D to -10.00D (0.50 steps after -8.00D) | 4 | Daily or up to 7 days extended | 13.3 | 6-pack blisters | 55 (methafilcon A) /6 |  |
| | | 8.3 8.9 | 14.2 | -0.25D to -8.00D (0.25 steps) | | | | | | |
| | | 8.8 | 14.3 | +0.25D to +8.00D (0.25 steps) | | | | | | |

Monthly Replacement Lenses

| | | | | | | | | | | |
|--|-----------------------|-----|------|--|-------------|---|------|--------------------|-------------------------------|--|
| Biofinity®  | comfilcon A/ 48% | 8.6 | 14.0 | +8.00D to -12.00D (0.50 steps after +/-6.00D) | SiHy (1) | Daily or 6 nights/ 7 days Extended | 160 | 6-pack blisters | SiH 48 (comfilcon A)/ 3 |  |
| Expressions® Accents  | methafilcon A/ 55% | 8.7 | 14.4 | +4.00D to -6.00D (0.25 steps) | 4 | Daily or up to 7 days extended | 37.6 | 6-pack blisters | No Trials | Colors available: Blue, Green, Violet, Aqua |
| Expressions® Colors  | methafilcon A/ 55% | 8.7 | 14.4 | +4.00D to -6.00D (0.25 steps) | 4 | Daily or up to 7 days extended | 37.6 | 6-pack blisters | No Trials | Colors available: Jade, Blue Topaz, Aqua, Blue, Green, Gray, Hazel, Brown |



PC TECHNOLOGY



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS














OPTIMIZED BALLAST DESIGN










UV BLOCKING

Sphere Lenses

| BRAND | MATERIAL/ H ₂ O CONTENT | BASE CURVE (mm) | DIA (mm) | SPHERE POWER | FDA Group | WEARING SCHEDULE | Dk/t (at -3.00D) | REVENUE CARTON SIZE | TRIALS LABELED AS/ PACK SIZE | FEATURES |
|--|---------------------------------------|--------------------|-------------|---|--------------|--------------------------------------|------------------------|---------------------------|------------------------------------|---|
| Frequency® 55  | Methafilcon A/ 55% | 8.4 | 14.2 | -0.25D to -8.00D (0.25 steps) | 4 | Daily or up to 7 days extended | 23.5 | 6-pack blisters | 55 (methafilcon A)/ 6 |   |
| | | 8.7 | 14.2 | +8.00D to -10.00D (0.50 steps after -8.00D) | | | | | | |
| | | 9.0 | 14.2 | -0.25D to -8.00D (0.25 steps) | | | | | | |
| Frequency® 55 Aspheric  | methafilcon A/ 55% | 8.7 | 14.4 | +8.00D to -10.00D (0.50 steps after -8.00D) | 4 | Daily or up to 7 days extended | 20.9 | 6-pack blisters | Aspheric (methafilcon A)/ 3 |    |
| | | 8.4 | 14.4 | -0.25D to -8.00D (0.25 steps) | | | | | | |
| Proclear® Sphere  | omafilcon A/ 62% | 8.6 | 14.2 | +20.00D to -20.00D (0.50 steps after +/-6.00D) | 2 | Daily | 52.3 | 6-pack | Proclear (omafilcon A)/ 3 |    |
| | | 8.2 | 14.2 | Plano to -10.00D (0.50 steps after -6.00D) | | | | | | |

Quarterly Replacement Lenses

| | | | | | | | | | | |
|--------------------------------------|-----------------------|-----------------------------|----------------------|--|---|--------------------------------------|--|-----------------------------|-----------|---|
| Edge® III ProActive | polymacon/ 38% | 8.4, 8.7 | 14.0 | -0.25D to -8.00D (0.50 steps after -6.00D) | 1 | Daily | 14 | single & 3 pk blister | No Trials |   |
| | | 8.7 | | +0.25D to +5.00D | | | | | | |
| Hydrasoft® Sphere Options | methafilcon B/ 55% | 8.6, 8.9 9.2 8.3, 8.6 | 15.0 15.0 14.2 | +10.00D to -20.00D (0.50 steps after -10.00D) | 4 | Daily or up to 7 days extended | 15.7 (15 Dia) 18.8 (14.2 Dia) | 4-pack vials | No Trials |  |
| Preference® | tetrafilcon A/ 43% | 8.4 8.7 | 14.4 | -0.25D to -10.00D (0.50 steps after -6.50D) | 1 | Daily or up to 7 days extended | 26.6 | 4-pack vials | No Trials |   |
| Preference® Standard | tetrafilcon A/ 43% | 8.3 8.6 | 14.0 14.0 | Plano to -10.00D (0.50 steps after -6.50D) | 1 | Daily | 13.3 | 4-pack vials | No Trials |   |
| | | 8.7 | 14.4 | +0.25D to -10.00D(0.25 steps) | | | | | | |



PC TECHNOLOGY



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS



OPTIMIZED BALLAST DESIGN



UV BLOCKING

Sphere Lenses

| BRAND | MATERIAL/ H ₂ O CONTENT | BASE CURVE (mm) | DIA (mm) | SPHERE POWER | FDA Group | WEARING SCHEDULE | Dk/t (at -3.00D) | REVENUE CARTON SIZE | TRIALS LABELED AS/ PACK SIZE | FEATURES |
|--|---------------------------------------|--|--------------------------------------|--|--------------|--------------------------------------|-----------------------------------|---------------------------|------------------------------------|--|
| Conventional Replacement Lenses | | | | | | | | | | |
| CooperClear® | tetrafilcon A/ 43% | 8.3 8.6 8.7 | 14.0 14.0 14.4 | +6.00D to -10.00D (0.50 steps after -6.50D) | 1 | Daily | 13.3 | single vial | No Trials | |
| CooperClear® FW | tetrafilcon A/ 43% | 8.4 8.7 | 14.4 | -0.25D to -10.00D (0.50 steps after -6.50D) | 1 | Daily or up to 7 days extended | 13.3 | single vial | No Trials | |
| Cooper® Prosthetic | polymacon/ 38% | 8.4 8.7 9.0 | 14.5 | +20.00D to -20.00D (0.50 steps after 10.00D) | 1 | Daily | | single vial | No Trials | Nine color combinations with open or closed pupils and black or clear backing |
| Edge® III ProActive XT | polymacon/ 38% | 8.4 8.7 9.0 | 14.0 | -0.25D to -8.00D (0.50 steps after -6.00D) | 1 | Daily | 14 | 2-pack blisters | No Trials | |
| Edge® III ProActive XT Plus | polymacon/ 38% | 8.7 | 14.0 | +0.25D to +5.00D | 1 | Daily | 14 | single vial | No Trials | |
| Hydrasoft® Standard | methafilcon B/ 55% | 8.6 8.9 9.2 8.3 8.6 | 15.0 15.0 15.0 14.2 14.2 | +10.00D to -20.00D (0.50 steps after -10.00D) | 4 | Daily or up to 7 days extended | 15.7 (15Dia) 18.8 (14.2) | single vial | No Trials | |
| Hydrasoft® XW Standard | methafilcon B/ 55% | 8.6 8.9 9.2 8.3 8.6 | 15.0 15.0 15.0 14.2 14.2 | +10.00D to -20.00D (0.50 steps after +/-10.00D) | 4 | Daily or up to 7 days extended | 31.3 | single vial | No Trials | |
| Hydrasoft® Custom Sphere | methafilcon B/ 55% | Call for availability including prism spheres and custom base curves and diameters | | | 4 | Daily | | single vial | No Trials | |



PC TECHNOLOGY



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS



OPTIMIZED BALLAST DESIGN



UV BLOCKING

Sphere Lenses

| BRAND | MATERIAL/ H ₂ O CONTENT | BASE CURVE (mm) | DIA (mm) | SPHERE POWER | FDA Group | WEARING SCHEDULE | Dk/t (at -3.00D) | REVENUE CARTON SIZE | TRIALS LABELED AS/ PACK SIZE | FEATURES |
|-----------------------------------|---------------------------------------|---------------------------------|--------------------------------------|--|-----------|--------------------------------|---------------------|---------------------|---------------------------------|---|
| Hydrasoft® Aphakic | methafilcon B/ 55% | 8.6 8.9 9.2 8.3 8.6 | 15.0 15.0 15.0 14.2 14.2 | +10.50D to +20.00D (0.50 steps) | 4 | Daily | | single vial | No Trials | |
| Hydrasoft® Aphakic Options | methafilcon B/ 55% | 8.6 8.9 9.2 8.3 8.6 | 15.0 15.0 15.0 14.2 14.2 | +10.50D to +20.00D (0.50 steps) | 4 | Daily or up to 7 days extended | | 4-pack vials | No Trials | |
| Hydrasoft® XW Aphakic | methafilcon B/ 55% | 8.6 8.9 9.2 8.3 8.6 | 15.0 15.0 15.0 14.2 14.2 | +10.50D to +20.00D (0.50 steps) | 4 | Daily or up to 7 days extended | | single vial | No Trials | |
| Hydron® ZERO 4® | polymacon/ 38% | 8.6 | 14.0 | -0.25D to -10.00D (0.50 steps after -6.00D) | 1 | Extended | 21 | single vial | No Trials | |
| Hydron® ZERO 6® | Polymacon/ 38% | 8.4 8.7 9.0 | 14.0 | +10.00D to -10.00D (0.50 steps after +5.00D and -6.00D) | 1 | Daily | 14 | single vial | No Trials | |
| Natural Touch® | polymacon/ 38% | 8.4 8.7 | 14.0 | Plano to -6.00D (0.25 steps) | 1 | Daily | 14 | single vial | No Trials | Colors available: Sophisticated Blue, Baby Blue, Hazel, Willow Green, Sultry Gray, Aqua Seas |
| Permalens® (Low Plus) | perfilcon A/ 71% | 8.0 8.3 8.6 | 14.0 | +0.25D to +8.00D (0.50 steps after +5.00D) | 4 | Extended | 12.5 | single vial | No Trials | |

PC
TECHNOLOGYVISIBILITY
TINTLOW
MODULUSNATURALLY
WETTABLEHIGH OXYGEN
TRANSMISSIBILITYASPHERIC
OPTICSOPTIMIZED
BALLAST DESIGNUV
BLOCKING

Sphere Lenses

| BRAND | MATERIAL/ H ₂ O CONTENT | BASE CURVE (mm) | DIA (mm) | SPHERE POWER | FDA Group | WEARING SCHEDULE | Dk/t (at -3.00D) | REVENUE CARTON SIZE | TRIALS LABELED AS/ PACK SIZE | FEATURES |
|-----------------------------------|---------------------------------------|-----------------------|-------------|---|--------------|--------------------------------------|------------------------|---------------------------|------------------------------------|--|
| Permalens® (Minus) | perfilcon A/ 71% | 7.7 | 13.5 | -0.25D to -20.00D (0.50 steps after -6.00D) | 4 | Extended | 34 | single vial | No Trials | |
| | | 8.0 | | | | | | | | |
| | | 8.3 | | | | | | | | |
| | | 8.6 | 14.2 | -0.25D to -10.00D (0.50 steps after -6.00D) | | | | | | |
| Permalens® Aphakic | perfilcon A/ 71% | 8.0 | 14.0 | +8.50D to +20.00D (0.50 steps) | 4 | Extended | 34 | single vial | No Trials | |
| | | 8.3 | 14.0 | | | | | | | |
| | | 8.3 | 14.5 | | | | | | | |
| | | 8.6 | 14.5 | | | | | | | |
| | | 8.9 | 14.5 | | | | | | | |
| Permalens® Therapeutic | perfilcon A/ 71% | 7.7 | 13.5 | Plano | 4 | Extended | 14.2 | single vial | No Trials | |
| | | 8.0 | 13.5 | | | | | | | |
| | | 8.3 | 13.5 | | | | | | | |
| | | 8.6 | 14.2 | | | | | | | |
| | | 9.0 | 15.0 | | | | | | | |
| Proclear® Sphere | omafilcon A/ 59% | 8.2 | 14.2 | +10.00D to -10.00D (0.50 steps after -8.50D and after +6.00D) | 2 | Daily | 47.1 | single vial | No Trials | |
| | | 8.5 | | | | | | | | |
| | | 8.8 | | | | | | | | |
| Vantage® | tetrafilcon A/ 43% | 8.3 | 14.0 | +6.00D to -10.00D (0.50 steps after -6.50D) | 1 | Daily | 13.3 | single vial | No Trials | |
| | | 8.6 | 14.0 | | | | | | | |
| | | 8.9 | 14.0 | | | | | | | |
| | | 8.7 | 14.4 | | | | | | | |
| | | 8.7 | 14.4 | | | | | | | |
| Vantage® Accents | tetrafilcon A/ 43% | 8.3 | 14.0 | Plano to -6.50D (0.25 steps) | 1 | Daily | 13.3 | single vial | No Trials | Enhancement colors available: Sky Blue, Violet Blue, Spring Green, Turquoise, Auburn, Misty Brown |
| | | 8.6 | 14.0 | | | | | | | |
| | | 8.7 | 14.4 | | | | | | | |
| | | 8.7 | 14.4 | | | | | | | |
| Vantage® Thin | tetrafilcon A/ 43% | 8.4 | 14.0 | Plano to -10.00D (0.50 steps after -6.50D) | 1 | Daily or up to 7 days extended | 26.6 | single vial | No Trials | |
| | | 8.7 | 14.4 | | | | | | | |
| Vantage® Thin Accents | tetrafilcon A/ 43% | 8.4 | 14.0 | Plano to -6.50D (0.25 steps) | 1 | Daily or up to 7 days extended | 26.6 | single vial | No Trials | Enhancement colors available: Sky Blue, Violet Blue, Spring Green, Turquoise, Auburn, Misty Brown |
| | | 8.7 | 14.4 | | | | | | | |



PC TECHNOLOGY



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS

























OPTIMIZED BALLAST DESIGN



UV BLOCKING

Toric Lenses

| BRAND | MATERIAL/ H ₂ O CONTENT | BASE CURVE (mm) | DIA (mm) | SPHERE POWER | CYLINDER POWER | AXIS | FDA Group | WEARING SCHEDULE | Dk/t (at -3.00D) | REVENUE CARTON SIZE | TRIALS LABELED AS/ PACK SIZE | FEATURES |
|---|---------------------------------------|----------------------|-------------|--|----------------------------------|----------------------------|-----------|--------------------------------|------------------|---------------------|-----------------------------------|---|
| 1 Day Disposable Lenses | | | | | | | | | | | | |
| ClearSight™ 1 Day Toric  | ocufilcon D/ 55% | 8.7 | 14.5 | Plano to -7.00D (0.50 steps after -6.00D) | -0.75 -1.25 | 180° 160° 90° 20° | 4 | Daily | 17.8 | 30-lens carton | 55% 1 Day Toric (ocufilconD)/ 5pk |    |
| 2- Week Replacement Lenses | | | | | | | | | | | | |
| Biomedics® Toric  | ocufilcon D/ 55% | 8.7 | 14.5 | +6.00D to -9.00D (0.50 steps after -6.00D and +5.00D) | -0.75 -1.25 -1.75 -2.25 | 10° to 180° (10° steps) | 4 | Daily or Extended | 17.8 | 6-pack blisters | Toric contact lens/ single |     |
| Vertex Toric®  | methafilcon A/ 55% | 8.6 | 14.4 | +6.00D to -8.00D (0.50 steps after -6.50D) | -0.75 -1.25 -1.75 -2.25 | 10° to 180° (10° steps) | 4 | Daily or up to 7 days extended | 17.9 | 6-pack blisters | Vertex toric/ single |   |
| Vertex Toric® XR  | methafilcon A/ 55% | 8.6 (labeled MED) | 14.4 | +6.00D to -8.00D (0.50 steps after -6.50D) | -2.75 -3.25 -3.75 | 5° to 180° (5° steps) | 4 | Daily or up to 7 days extended | 17.9 | 6-pack blisters | Toric 55/ single |   |
| Monthly Replacement Lenses | | | | | | | | | | | | |
| Biofinity® Toric  | comfilcon A/ 48% | 8.7 | 14.5 | +8.00 to -10.00D (0.50 steps after +/-6.00) | -0.75 -1.25 -1.75 -2.25 | 10° to 180° (10° steps) | SiHy (1) | Daily or up to 7 days extended | 116 | 6-pack blisters | Comfilcon A Dk 128 Toric/ single |      |
| Frequency® 55 Toric  | methafilcon B/ 55% | 8.4 8.7 | 14.4 | +6.00D to -8.00D (0.50 steps after -6.50D) | -0.75 -1.25 -1.75 -2.25 | 10° to 180° (10° steps) | 4 | Daily or up to 7 days extended | 15.7 | 6-pack blisters | Frequency 55 toric/ single |   |
| Frequency® 55 Toric XR  | methafilcon B/ 55% | 8.4 8.7 | 14.4 | +6.00D to -8.00D (0.50 steps after -6.50D) | -2.75 to -5.75 (0.50D steps) | 5° to 180° (5° steps) | 4 | Daily or up to 7 days extended | 15.7 | 6-pack blisters | Frequency 55 toric/ single |   |



PC TECHNOLOGY



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS


















OPTIMIZED BALLAST DESIGN



UV BLOCKING

Toric Lenses

| BRAND | MATERIAL/ H ₂ O CONTENT | BASE CURVE (mm) | DIA (mm) | SPHERE POWER | CYLINDER POWER | AXIS | FDA Group | WEARING SCHEDULE | Dk/t (at -3.00D) | REVENUE CARTON SIZE | TRIALS LABELED AS/ PACK SIZE | FEATURES |
|---|---------------------------------------|--------------------------|------------------------------|---|--|----------------------------|-----------|--------------------------------|---------------------|---------------------|---------------------------------|---|
| Proclear® Toric  | omafilcon A/ 62% | 8.8 8.4 | 14.4 | +6.00D to -8.00D (0.50 steps after -6.50D) | -0.75 -1.25 -1.75 -2.25 | 10° to 180° (10° steps) | 2 | Daily | 30.9 | 6-pack blisters | Proclear toric/ single |    |
| Proclear® Toric XR  | omafilcon A/ 59% | 8.8 8.4 | 14.4 | +6.25D to +10.00D (0.50 steps after +6.50D) -8.50D to -10.00D (0.50 steps) | -0.75 -1.25 -1.75 -2.25 | 5° to 180° (5° steps) | 2 | Daily | 30.9 | 6-pack blisters | Proclear toric/ single |    |
| | | 8.8 8.4 | 14.4 | +10.00D to -10.00D (0.50 steps after +/-6.50D) | -2.75 to -5.75 (0.50D steps) | | | | | | | |
| Quarterly Replacement Lenses | | | | | | | | | | | | |
| Cooper Toric™ | tetrafilcon A/ 43% | 8.4 8.7 | 14.4 | +6.00D to -8.00D (0.50 steps) | -0.75,-1.25 -1.75, -2.25 | 10° to 180° (10° steps) | 1 | Daily | 10.3 | Single vial | No Trials |   |
| Hydrasoft® Toric Options | methafilcon B/ 55% | 8.3 8.6 8.9 9.2 | 15.0 15.0 15.0 15.0 | +20.00D to -20.00D (0.25 steps) | -0.50 to -10.00 (-0.50 cylinder not available on sphere powers over +10.00) | any axis | 4 | Daily or up to 7 days extended | 14.5 | 4-pack vials | No Trials |  |
| | | 8.3 8.6 | 14.2 14.2 | +10.00D to -20.00D (0.25 steps) | -0.50 to -10.00 (0.25steps) | | | | | | | |
| Preference® Toric | tetrafilcon A/ 43% | 8.4 8.7 | 14.4 | +6.00D to -8.00D (0.50 steps after -6.50D) | -0.75, -1.25 -1.75, -2.25 | 5° to 180° (5° steps) | 1 | Daily | 10.3 | 4-pack vials | No Trials |   |
| Preference® Toric XR | tetrafilcon A/ 43% | 8.4 8.7 | 14.4 | +6.00D to -9.50D (0.50 steps after -6.50D) <small>*For sphere powers unique to XR cylinder powers below -2.75 are available</small> | -2.75 to -9.75 (0.25 steps) | 5° to 180° (5° steps) | 1 | Daily | 10.3 | 4-pack vials | Preference toric/ single |   |

PC
TECHNOLOGYVISIBILITY
TINTLOW
MODULUSNATURALLY
WETTABLEHIGH OXYGEN
TRANSMISSIBILITYASPHERIC
OPTICSOPTIMIZED
BALLAST DESIGNUV
BLOCKING

Toric Lenses

| BRAND | MATERIAL/ H ₂ O CONTENT | BASE CURVE (mm) | DIA (mm) | SPHERE POWER | CYLINDER POWER | AXIS | FDA Group | WEARING SCHEDULE | Dk/t (at -3.00D) | REVENUE CARTON SIZE | TRIALS LABELED AS/ PACK SIZE | FEATURES |
|---|---------------------------------------|---|--|--------------------|--------------------|-----------------------|--------------|--------------------------------------|------------------------|---------------------------|------------------------------------|----------|
| Conventional Replacement Lenses | | | | | | | | | | | | |
| Hydrasoft® Toric DW/XW Division I Standard | methafilcon B/ 55% | 8.9 8.6 | 15.0 14.2 | +3.00D to -6.00D | -0.75 to -2.00 | 180 ± 20° 90 ± 20° | 4 | Daily or up to 7 days extended | 14.5 | single vial | No Trials | |
| Hydrasoft® Toric DW/XW Division II Specialty | methafilcon B/ 55% | 8.6 8.9 8.3 8.6 | 15.0 15.0 14.2 14.2 | +10.00D to -20.00D | -0.75 to -5.00 | any axis | 4 | Daily or up to 7 days extended | 14.5 | single vial | No Trials | |
| Hydrasoft® Toric DW/XW Division III Custom | methafilcon B/ 55% | 8.3 8.6 8.9 9.2 8.3 8.6 | 15.0 15.0 15.0 15.0 14.2 14.2 | +10.00D to -20.00D | -0.50 to -10.00 | any axis | 4 | Daily or up to 7 days extended | 14.5 | single vial | No Trials | |
| Hydrasoft® Toric DW/XW Division III Aphakic | methafilcon B/ 55% | 8.6 8.9 | 15.0 | +10.25D to +20.00D | -0.75 to -10.00 | any axis | 4 | Daily or up to 7 days extended | 14.5 | single vial | No Trials | |
| Hydrasoft® Toric DW/XW Division IV Custom | methafilcon B/ 55% | Call for availability on special parameters | | | | | 4 | Daily or up to 7 days extended | 14.5 | single vial | No Trials | |



PC TECHNOLOGY



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS












OPTIMIZED BALLAST DESIGN



UV BLOCKING

Multifocal Lenses

| BRAND | MATERIAL/ H ₂ O CONTENT | BASE CURVE (mm) | DIA (mm) | SPHERE POWER | ADD POWER | DESIGN | FDA Group | WEARING SCHEDULE | Dk/t (at -3.00) | REVENUE CARTON SIZE | TRIALS LABELED AS/ PACK SIZE | FEATURES |
|---|---------------------------------------|--------------------|-------------|--|---|---|-----------|--------------------------------|--------------------|---------------------|--------------------------------------|---|
| 2- Week Replacement Lenses | | | | | | | | | | | | |
|  ProcLEAR® EP | omafilcon A/ 60% | 8.7 | 14.4 | +6.00D to -8.00D (0.50 steps after -6.50D) | LOW ADD (For patients with up to +1.25 ADD) | Center spherical distance, progressive intermediate and near aspheric | 2 | Daily | 33 | 6-pack blisters | 60% EP/ 3 |  |
| Monthly Replacement Lenses | | | | | | | | | | | | |
|  Frequency® 55 Multifocal | methafilcon A/ 55% | 8.7 | 14.4 | +4.00D to -6.00D (0.25 steps) | +1.00, +1.50, +2.00, +2.50 | D lens N lens | 4 | Daily or up to 7 days extended | 23.5 | 6-pack blisters | Frequency 55 Multifocal/ single |  |
|  ProcLEAR® Multifocal | omafilcon A/ 62% | 8.7 | 14.4 | +6.00D to -8.00D (0.50 steps after -6.50D) | +1.00, +1.50, +2.00, +2.50 | D lens N lens | 2 | Daily | 21.3 | 6-pack blisters | ProcLEAR Multifocal/ single |  |
|  ProcLEAR® Multifocal XR | omafilcon A/ 59% | 8.7 | 14.4 | +6.25D to +20.00D -8.50D to -20.00D (0.50 steps after +/-6.50D) | +1.00, +1.50, +2.00, +2.50 | D lens N lens | 2 | Daily | 21.3 | 6-pack blisters | ProcLEAR Multifocal/ single |  |
| | | 8.7 | 14.4 | +20.00D to -20.00D (0.50 steps after +/- 6.50D) | +3.00, +3.50, +4.00 | | | | | | | |
| | | 8.4 | 14.4 | +20.00D to -20.00D (0.50 steps after +/- 6.50D) | +1.00, +1.50,+2.00, +2.50,+3.00, +3.50, +4.00 | | | | | | | |
|  ProcLEAR® Multifocal Toric | Omafilcon A/ 59% | 8.8 8.4 | 14.4 | +20.00D to -20.00D (0.50 steps after +/- 6.50D) CLY: -0.75 to -5.75 (0.50 steps) AXIS: 5° to 180° (5° steps) | +1.00 to +4.00 (0.50 steps) | D lens N lens | 2 | Daily | 48.6 | 6-pack blisters | ProcLEAR Multifocal Toric/ single |  |

PC
TECHNOLOGYVISIBILITY
TINTLOW
MODULUSNATURALLY
WETTABLEHIGH OXYGEN
TRANSMISSIBILITYASPHERIC
OPTICSOPTIMIZED
BALLAST DESIGNUV
BLOCKING**“IT’S OKAY” GUARANTEE:**

CooperVision's "It's Okay" Guarantee allows exchanges or returns of all the CooperVision lenses listed below during the first 100 days following the original invoice. While products NOT listed below cannot be returned under the "It's Okay" guarantee, they may still be able to be returned under a different CooperVision return policy.

Product is invoiced at the full price and full credit is issued for all products returned within 100 days of the original invoice. All unopened product being returned for credit must be intact and in re-saleable condition; free of writing, stickers, or markings of any kind. Any product not received in re-saleable condition will not be credited and will be destroyed. Patient cancellation credit made within 100 days of original lens purchase will be 100% of the original lens price (as long as they are returned in re-saleable condition).

For full credit on opened products, 90-pack returns/exchanges must contain at least 25 unopened lenses in intact strips, 6-pack returns/exchanges must contain at least 5 unopened lenses, and 4-pack returns/exchanges must contain at least 3 unopened

ELIGIBLE BRANDS:

Avaira®
Biofinity®
Biomedics®
ClearSight™
Hydrasoft®
Proclear®
Vertex®

MADE-TO-ORDER LENSES:

CooperVision offers some lenses that are made-to-order and will ship within 10 days from the time of order. These lenses are Proclear Toric XR, Proclear Multifocal XR, Proclear Multifocal Toric, Vertex Toric XR, Frequency 55 Toric XR, and Preference Toric XR. Hydrasoft Toric and/or Options will ship within 14-21 days from the time of order.

UV NOTICE:

UV-absorbing contact lenses are NOT substitutes for protective UV-absorbing eyewear, such as UV-absorbing goggles or sunglasses, as they do not completely cover the eye and surrounding area. Patients should continue to use UV-absorbing eyewear as directed.